



INTEGRATED
EMBOLIC
PROTECTION

VANGUARD IEP[®]

EMBOLIC PROTECTION WITH EASE.

EVERY TIME.




Vanguard IEP Peripheral Balloon Angioplasty System with Integrated Embolic Protection

The **ONLY** lower extremity balloon angioplasty system that includes an integrated filter on the same catheter

- The **ONLY** filter adjustable to patient anatomy
- One device equals less exchanges
- Sheathless filter design


INTEGRATED EMBOLIC PROTECTION™

Vanguard IEP



The Vanguard IEP Peripheral Balloon Angioplasty System with Integrated Embolic Protection

The Vanguard IEP System combines an over-the-wire, percutaneous transluminal angioplasty (PTA) catheter with a nitinol embolic protection filter system, positioned distally to the angioplasty balloon.

- Single step filter deployment
- Minimal balloon-to-filter distance (5mm)
- .018" wire compatible

Available Size Matrix

VANGUARD IEP System Sizes		Available Catheter Length (cm)	Balloon Length (mm) Reference Code			Sheath Compatibility	Filter Length (mm)	Filter Diameter (mm)
			80	120	200			
Balloon OD (mm)	5.0	135	V1-5080-135-01	V1-5120-135-01	V1-5200-135-01	6F	16	Adjustable up to 8 mm
	6.0	135	V1-6080-135-01	V1-6120-135-01	V1-6200-135-01	6F		

The Vanguard IEP Peripheral Balloon Angioplasty System with Integrated Embolic Protection is indicated for peripheral vascular percutaneous transluminal angioplasty (PTA) and capture and removal of embolic material during angioplasty for the femoral, iliac, popliteal and profunda arteries.

The Vanguard IEP is CE marked and is available for sale in select European markets. Vanguard IEP, Contego Medical, IEP and Integrated Embolic Protection are trademarks of Contego Medical. © 2017 Contego Medical, LLC. All rights reserved.

Contego Medical, LLC

3921 Sunset Ridge Road, Suite 102, Raleigh, NC 27607 USA | PH 01 919 459 7250
 HealthLink Europe BV 5215 MC's Hertogenbosch, The Netherlands | PH 0031 733032517
www.contegomedical.com MKT-1084-04 DCO#0407 March 2018

